

Hi Skin

Skin Care Platform

Addresses All Skincare Needs

Multifunctional Skin Care Platform 9 in 1




FINE LINES + WRINKLES




ELASTICITY + FIRMNESS



EVEN TONE + VIBRANCY



SKIN TEXTURE



BROWN SPOTS



OILY + CONGESTED SKIN



ENLARGED PORES



NINE DIFFERENT HANDLES



Hydro Dermabrasion



Polymer Atomizing Spray



Ultrasound Import



Bipolar RF

Multi-pole Microcurrent



Bipolar Microcurrent



Cold Hammer



Ultrasonic Vibration



Mousse Bubbles




DEEP CLEANING




Hydro Dermabrasion


Open pores, soften the stratum corneum, deeply cleanse the skin, deeply moisturize, remove blackheads, improve the dark yellow skin, and make the skin clean, white and supple.



Solution is sprayed through the inlet of nozzle tip. It rotates clockwise.



The waste and remaining solution are sucked to container through inlets.



Solution are supplied to skin surface through nozzle tip.



Hard turbine hydro cleaning head



Microdermabrasion head



Soft silicone hydro cleaning head



Pipe cleaning head



DEEP CLEANING

Mousse Bubbles Handle

The handle produces dense foam containing a large amount of active oxygen components. Using the high oxidative properties of active oxygen, the hydroxide ions can quickly penetrate into the skin, so as to achieve the effect of cleaning, moisturizing and rejuvenating the skin.



Ultrasonic Vibration Handle

Remove dead skin by high-frequency vibration waves, remove blackheads, soften skin keratinocytes, and deeply cleanse skin dirt.




MOISTURIZING



Ultrasound Import

At a frequency of 1 million to 3 million vibrations per second to quickly open the pores, allowing the essence to penetrate into the basal layer of the skin



Polymer Atomizing Spray

Combined with liquid products, the liquid is converted into mist-like small molecules by high-pressure atomization, and then quickly poured into the dermis by high-pressure spraying to achieve a deep hydrating effect



TIGHTENING & LIFTING



Microcurrent Handles

Microcurrent uses low-voltage electricity to stimulate muscle, adenosine triphosphate (ATP) cell growth, and collagen development in the dermis on the face.

Biopolar RF Handle

RF treatment is a non-invasive way to achieve skin tightening and make you look younger. The RF energy will heat the deep skin layers (dermis) of your body to produce collagen.





Cold & Warm Hammer Handle

Cold Hammer 0°C: Ice skin, shrink pores, calm skin

Warm Hammer 45°C: increase the temperature of the skin surface, accelerate blood circulation, and increase the absorption rate when used with products.



SCOPE OF TREATMENT

Deep Cleaning

Remove Blackheads, Mites Clear, Remove grease dust,
Remove Dead Skin

Skin Rejuvenation

Skin Moisturizing, Shrink Pores, Essence Import, Skin
Whitening, Acne Treatment

Skin Tightening

Face Lifting, Eye Lifting, Wrinkle removal





SPECIFICATION

Power	600w
RF frequency	3Mhz
Screen	8.4inch
Vacuum	0-90kp
Flow Capacity	30L Max
Probe of dermabrasion	11pics
Voltage	100V-240V 50Hz-60Hz